



# RelaxView Healthcare Lighting System

## RelaxView Ambient Ceiling Illuminator RVAC 1x1 Flange (Surface Mount)

RelaxView Ambient Ceiling Illuminator RVAC 1x1 Flange (Surface Mount)

The Day-Brite / CFI RelaxView Healthcare Lighting Solution provides a complete digital lighting solution benefiting both healthcare professional and patient. The LED Graphic Panels may be mounted on the ceiling and wall as a distraction during long procedures. While healthcare professionals have general ambient lighting to aid in the performance of their tasks. Wall dimmers and power distribution modules provide additional control of the environment for the task at hand as well as navigation throughout the room.

### Product data

General information	
Lamp Source Code	LED
Colour Temperature	3500 4200 K K
CRI	0
Photometric Distribution Type	Direct
Mounting	Surface T-Bar / T-Grid
Mechanical and housing	
Fixture Size	1 ft x 1 ft
Shape	Square
Initial performance (IEC compliant)	
Luminaire Lumens	0

Product data	
Order product name	LED Ambient Ceiling Illuminator RVAC
Order code	RVAC11FS
Material Nr. (12NC)	DAY70000373
Catalog Number	RVAC11FS
Catalog Number Description	RelaxView Ambient Ceiling Illuminator RVAC 1x1 Flange (Surface Mount)

## RelaxView Healthcare Lighting System

