



# CareWell Bed Light

## Assisted Living Patient Bed Light

### Assisted Living Patient Bed Light

CareWell is a value oriented solution for resident apartments within assisted living facilities. The varied profiles and finishes of the fast-growth poplar wood shroud provide a wide range of styles that provide a residential feel. Dual lamp compartments deliver ambient and reading illumination while eliminating debilitating direct glare. CareWell is easy to install and maintain.

#### Product data

General Information	
Lamp Source Code	FLU
Lamp Watts	32
Lamp Envelope Code	T-8
Number of Lamps	2
Photometric Distribution Type	General Diffuse
Operating and Electrical	
Voltages	120/277/347 V
Mechanical and Housing	
Fixture Size	48 in

Shape	Linear
Product Data	
Order product name	Assisted Living Patient Bed Light
Commercial Code	ALB4T8-SYM-11-120-SWH
Order code	ALB4T8-SYM-11-120-SWH
Catalog Number Description	Assisted Living Patient Bed Light
Service / Spec Smart	-
Material Nr. (12NC)	ALK00004021